

What is Derma Roller and Meso Roller?

Derma roller, which is also be known as meso roller, is an outstanding cosmetic product of the modern technology, it consists of a roller containing 540 micro-needles and a handle; applies to the skin and stimulates the body into producing collagen and elastin as a repair mechanism. This large number of micro-needles enables less roll but open more pores for the treatment serum to penetrate more into the skin, definitely helpful in stretch marks, scars, wrinkles, anti-aging, smoker's skin, skin rejuvenation and glowing, etc.

A cosmetic handheld device which consists of a roller containing 192 micro-needles and a handle which is used to treat a range of skin conditions. It is applied to the skin where it stimulates the body into producing collagen and elastin as a repair mechanism. (Our DermaRoller has 540 micro-needles)

Derma Roller and Meso Roller are the same thing, just in different terms. Our Derma Roller contains 540 micro-needles, enabling to roll less but open more pores for the treatment serum to penetrate more into the skin. It helps in stretch marks, scars, wrinkles, anti-aging, smoker's skin, skin rejuvenation and glowing.

Use Derma Roller or Meso Roller for :

- Stretch Mark Removal
- Acne / Scar / De-Pigmentation
- Anti Wrinkle Treatment
- Anti Aging Treatment
- Hair Loss Treatment or Hair Restoration
- Scar Removal
- Cellulite Treatment and Cellulite Reduction
- Smoker's skin/ Sun Damage Skin/ Uneven Complexion
- Glowing/ Balancing Treatment

How do you use a Derma Roller?

You can use a derma roller at the comfort of your own home. Please visit links below for instructions. Derma Roller can be used on any part of your body, but we do not recommend using it on eyelid.

How long will I see results and is it safe to use a Derma Roller?

The result depends on the condition of your skin and how often you do the treatment on the area. Results will be seen within a week for face if you do at least 3 – 5 treatment a week. Stretch mark and cellulite treatment will take about 2 weeks if you perform treatment 3 – 5 times a week. It is very safe to use a Derma Roller. It is FDA approved and it is manufactured and quality control by FDA standards.

Does it hurt using a Derma Roller?

No, you will feel slightly tingling on the face, but not a pain on your body. Light pressure rolling on the skin and none of this should draw blood.



Recommendation on Use of Derma Roller

For light Stretch Mark removal, we recommend using 1mm derma roller 5 times/week.
For deep Stretch Mark removal, we recommend using 1.5mm 3 times/week.

For Cellulite removal, we recommend using 1mm derma roller 5 times/week. Or 1.5mm derma roller 5 times/week.

For Anti-Aging, Wrinkle, Acne/Scars/De-Pigmentation, Glowing, Anti-Sagging, we recommend using 0.5mm 3 to 5 times/week.